

NANOCOLOR® CSB 15000 / COD / DCO / DQO



Method(e):


0281
1.0–15.0 g/L O₂

20–80 %:


1.0 3.8 12.2 15.0 g/L O₂




REF 985 028




Probe / Sample
Échant. / Muestra




Schütteln / Shake
Agiter / Agitar




160 °C / 30 min
148 °C / 2 h



Schütteln / Shake
Agiter / Agitar



Abkühlen / Cool
Refroidir / Enfriar



Säubern / Clean
Nettoyer / Limpiar



Messen / Measure
Mesurer / Medir