

visicolor® ECO Alkalinität TA

Alkalinity TA / Alcalinité TA / Alkalinidad TA


Test 5-04

REF 931 204


436/585 nm

Method(e) / Método


5041	0.3–14 °d
5042	0.10–5.00 mmol/L H ⁺
5043	0.4–17.5 °e
5044	0.5–25.0 °f
5045	5–250 mg/L CaCO ₃




Probe / sample
échantillon / muestra




1
Küvetten spülen
Rinse test tubes
Rincer les éprouvettes
Lavar los tubos



2
5 mL Probe
Sample
Echantillon
Muestra




3
200 µL TA-1




4
1 NANOFIX TA-2



5
Schütteln
Shake
Agiter
Agitar



6
Säubern
Clean
Nettoyer
Limpiar




7
2'00 min



8
Messung
Measurement
Mesure
Medición

Umrechnungstabelle / Conversion table / Tableau de conversion / Tabla de conversión

°d	°e	°f	mg/L CaCO ₃	mmol/L H ⁺	gpg
1	1.3	1.8	18	0.36	1
2	2.5	3.6	36	0.72	2
3	3.8	5.4	54	1.08	3
4	5.0	7.1	71	1.42	4
5	6.3	8.9	89	1.78	5
6	7.5	10.7	107	2.14	6
7	8.8	12.5	125	2.50	7
8	10.0	14.3	143	2.86	8
9	11.3	16.1	161	3.22	9
10	12.5	17.8	178	3.56	10



Meerwasser
Sea water
Eau de mer
Agua de mar