

CSB 1500 Hg-frei

COD Hg-free/DCO sans Hg/DQO sin Hg

Test 0-29


REF 963 029




Rundküvettentest/ Tube test

Test en cuves rondes/ Prueba en cubeta cilíndrica


100–1500 mg/L




2.0 mL Probe
Sample
Echantillon
Muestra




Sicherheitsgefäß
Safety bottle
Récipient de sécurité
Recipiente de seguridad




148 °C/2 h




Schütteln
Shake
Agiter
Agitar



Abkühlen lassen
Allow to cool
Laisser refroidir
Dejar enfriar



Säubern
Clean
Nettoyer
Limpiar



Messung
Measurement
Mesure
Medición